

# FREE PROSTATE SPECIFIC ANTIGEN (F-PSA) ENZYME IMMUNOASSAY TEST KIT

## Catalog Number: BC-1021



**BioCheck, Inc**  
323 Vintage Park Dr.  
Foster City, CA 94404

### Enzyme Immunoassay for the Quantitative Determination of Free Prostate Specific Antigen (f-PSA) in Human Serum

#### FOR RESEARCH USE ONLY

Store at 2 to 8°C.

#### PROPRIETARY AND COMMON NAMES

BioCheck f-PSA Enzyme Immunoassay

#### INTENDED USE

For the quantitative determination of the Cancer Antigen f-PSA concentration in human serum.

#### INTRODUCTION

Human Prostate Specific Antigen (PSA) is a 33 kD serine proteinase which, in human serum, is predominantly bound to alpha 1-antichymotrypsin (PSA-ACT) and alpha 2-macroglobulin (PSA-AMG). Trace amounts of alpha 1-antitrypsin and inter-alpha trypsin inhibitor bound to PSA can also be found. Any remaining PSA is in the free form (f-PSA).<sup>1,3</sup> Current methods of screening men for prostate cancer utilize the detection of the major PSA-ACT form. Levels of 4.0 ng/ml or higher are strong indicators of the possibility of prostatic cancer.<sup>4</sup> However, elevated serum PSA levels have also been attributed to benign prostatic hyperplasia and prostatitis, leading to a large percentage of false positive screening results.<sup>5</sup> A potential solution to this problem involves the determination of free PSA levels.<sup>6-17</sup> Preliminary studies have suggested that the percentage of free PSA is lower in patients with prostate cancer than those with benign prostatic hyperplasia.<sup>2</sup> Thus, the measurement of free serum PSA in conjunction with total PSA, can improve specificity of prostate cancer screening in selected men with elevated total serum PSA levels, which would subsequently reduce unnecessary prostate biopsies with minimal effects on cancer detection rates.<sup>6</sup>

#### PRINCIPLE OF THE TEST

The f-PSA ELISA test is a solid phase two-site immunoassay. An anti-f-PSA monoclonal antibody is coated on the surface of the microtiter wells and a rabbit anti-PSA antibody labeled with horseradish peroxidase is used as the tracer. The f-PSA molecules present in the standard solution or sera are "sandwiched" between the two antibodies. Following the formation of the coated antibody-antigen- antibody-enzyme complex, the unbound antibody-enzyme tracers are removed by washing. The horseradish peroxidase activity bound in the wells is then assayed by a colorimetric reaction. The intensity of the color formed is proportional to the concentration of f-PSA present in the sample.

#### REAGENTS

##### Materials provided with the kits:

- Murine Monoclonal Anti-Free PSA coated microtiter plate with 96 wells.

- Reference standards containing 0, 1.0, 2.5, 5.0, 10.0, and 25.0 ng/ml f-PSA, 1 ml/vial, lyophilized. 1 set.
- f-PSA Zero Buffer, 13 ml.
- Enzyme Conjugate Reagent, 22 ml.
- TMB Reagent (One-Step), 11 ml.
- Stop Solution (1N HCl), 11 ml.

##### Materials required but not provided:

- Precision pipettes: 100  $\mu$ l, 200  $\mu$ l, and 1.0 ml.
- Disposable pipette tips.
- Distilled water.
- Vortex mixer or equivalent.
- Absorbent paper or paper towels.
- Microtiter plate shaker
- Graph paper.
- A microtiter plate reader with a bandwidth of 10nm or less and an optical density range of 0-2 OD or greater at 450 nm

#### SPECIMEN COLLECTION AND PREPARATION

Serum should be prepared from a whole blood specimen obtained by acceptable medical techniques. This kit is for use with serum samples without additives only.

#### STORAGE OF TEST KIT AND INSTRUMENTATION

Unopened test kits should be stored at 2-8°C upon receipt and the microtiter plate should be kept in a sealed bag with desiccants to minimize exposure to damp air. Opened test kits will remain stable until the expiration date shown, provided it is stored as described above. A microtiter plate reader with a bandwidth of 10 nm or less and an optical density range of 0-2 OD or greater at 450 nm wavelength is acceptable for use in absorbance measurement.

#### REAGENT PREPARATION

1. All reagents should be allowed to reach room temperature (18-25°C) before use.
2. Reconstitute each lyophilized standard with 1.0 ml distilled water. Allow the reconstituted material to stand for at least 20 minutes and mix gently. Reconstituted standards will be stable for up to 30 days when stored sealed at 2-8°C.

#### ASSAY PROCEDURE

1. Secure the desired number of coated wells in the holder.
2. Dispense 100  $\mu$ l of standards, specimens, and controls into appropriate wells.
3. Dispense 100  $\mu$ l of f-PSA Zero Buffer into each well.
4. Gently mix for 30 seconds.
5. Incubate at room temperature (18-25°C) for 60 minutes.
6. Remove the incubation mixture by emptying plate contents into a suitable waste container.
7. Rinse and empty the microtiter wells 5 times with distilled or deionized water. (Please do not use tap water.)
8. Strike the wells sharply onto absorbent paper or paper towels to remove all residual water droplets.
9. Dispense 200  $\mu$ l of Enzyme Conjugate Reagent into each well. Gently mix for 10 seconds.
10. Incubate at room temperature for 60 minutes.
11. Remove the incubation mixture by emptying plate contents into a suitable waste container.

12. Rinse and empty the microtiter wells 5 times with distilled or deionized water. (Please do not use tap water.)
13. Strike the wells sharply onto absorbent paper to remove residual water droplets.
14. Dispense 100  $\mu$ l of TMB Reagent into each well. Gently mix for 10 seconds.
15. Incubate at room temperature for 20 minutes.
16. Stop the reaction by adding 100  $\mu$ l of Stop Solution to each well.
17. Gently mix for 30 seconds. ***It is important to make sure that all the blue color changes to yellow color completely.***
18. Using a microtiter plate reader, read the optical density at 450 nm ***within 15 minutes.***

**Important Notes:**

1. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
2. It is recommended that if manual pipetting is used, no more than 32 wells be used for each assay run, since pipetting of all standards, specimens and controls should be completed within 3 minutes. A full plate of 96 wells may be used if automated pipetting is available.
3. Duplication of all standards and specimens, although not required, is recommended.

**CALCULATION OF RESULTS**

1. Calculate the average absorbance values ( $A_{450}$ ) for each set of reference standards, control, and samples.
2. Construct a standard curve by plotting the mean absorbance obtained for each reference standard against its concentration in ng/ml on linear graph paper, with absorbance on the vertical (y) axis and concentration on the horizontal (x) axis.
3. Using the mean absorbance value for each sample, determine the corresponding concentration of f-PSA in ng/ml from the standard curve.


**EXAMPLE OF STANDARD CURVE**

Results of a typical standard run with optical density readings at 450nm shown in the Y axis against f-PSA concentrations shown in the X axis. This standard curve is for the purpose of illustration only, and should not be used to calculate unknowns. Each user should obtain his or her own data and standard curve.

f-PSA (ng/ml)	Absorbance (450 nm)
0	0.058
1.0	0.222
2.5	0.447
5.0	0.800
10.0	1.400
25.0	2.817

**EXPECTED VALUES**

As discussed in the introduction, the important diagnostic parameter



is not the level of free PSA, but rather the ratio of free PSA to total PSA. Percent free-PSA offered the greatest advantage to the total PSA test when the total PSA values were between 3.0 and 10.0 ng/ml.<sup>14</sup>

***For a given patient sample, different commercial test kits of total PSA and free-PSA may give different values of total PSA and free-PSA. Users should keep this in mind while calculating the percentage. The following information is cited from References 6,7,10,11,13,14-17.***

For total PSA levels between 3.0 and 4.0 ng/ml, using a 19% cutoff point for percent free-PSA would result in detection of 90% of all cancers.<sup>14</sup> For total PSA levels between 4.1 and 10.0 ng/ml, the most appropriate cutoff point for free-PSA is 24%. At this cutoff point, 95% of the cancers would be detected.<sup>14</sup>

With respect to free PSA levels and prostate volume, the available information is again limited. Catalona et al were the first to demonstrate the importance of prostate size in selecting the cutoff value for percent free-PSA.<sup>13</sup> In their study, men with prostate cancer and a prostate volume of 400 cc or less had a median free-to-total PSA proportion of 0.092 (9.2%), a value statistically lower than the 0.159 (15.9%) found for patients with prostate cancer and a gland >40.0 cc. Yemoto et al, in a recent study of 200 men, showed no correlation between percent free-PSA and prostate volume.<sup>16</sup>

Data from several studies have demonstrated an inverse relationship between percent free-PSA and total PSA. This observation suggests that higher PSA levels are more commonly associated with lower percent free-PSA values and these men most frequently have more aggressive or advanced prostate cancers.<sup>2,14,17</sup>

**SENSITIVITY**

The sensitivity of this kit is estimated to be 0.1 ng/ml.

**LIMITATIONS OF THE PROCEDURE**

1. Reliable and reproducible results will be obtained when the assay procedure is carried out with a complete understanding of the package insert instructions and with adherence to good laboratory practice.

2. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
3. Serum samples demonstrating gross lipemia, gross hemolysis, or turbidity should not be used with this test.
4. The results obtained from the use of this kit should be used only as an adjunct to other diagnostic procedures and information available to the physician.

## REFERENCES

- 1 Christensson A, Bjork T, Nilsson O, et al. Serum Prostate Specific Antigen Complexed to  $\alpha$ 1-Antichymotrypsin As An Indicator of Prostate Cancer. *J. of Urol.* 150:100-105; 1993.
- 2 Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha$ -1-antichymotrypsin. *Clin Chem* 37:1618-1625, 1991.
- 3 Stenman U-H, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O. A complex between prostate-specific antigen and  $\alpha$ -1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: assay of the complex improves clinical sensitivity for cancer. *Cancer Res.* 51:222-226, 1991.
- 4 Catalona WJ, Smith DS, Ratliff TL, Basler JW. Detection of organ-confined prostate cancer is increased through prostate-specific antigen-based screening. *JAMA.* 270:948-954, 1993.
- 5 Stamey TA, Yang N Hay AR, McNeal JE, Freiha, FS, Redwine E. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. *N Engl J Med.* 317-:909-916, 1987.
- 6 Junker R, Brandt B, Zechel C, and Assmann, G. Comparison of Prostate-specific antigen (PSA) measured by four combinations of free-PSA and total PSA assays. *Clinical Chemistry* 43:1588-1594, 1997.
- 7 Luderer AA, Chen Y-T, Soriano TF, et al. Measurements of the proportion of free to total prostate-specific antigen improves diagnostic performance of prostate-specific antigen in the diagnostic gray zone of total prostate-specific antigen. *Urology* 46:187-194, 1995.
- 8 Carter HB, Pearson JD, Metter J, et al. Longitudinal evaluation of prostate-specific antigen levels in men with and without prostate disease. *JAMA.* 267:2215-2220, 1992.
- 9 Smith DS, Catalona WJ. Rate of change in serum prostate-specific antigen levels as a method for prostate cancer detection. *J Urol.* 152:1163-1167, 1994.
- 10 Benson MC, Whang IS, Pantuck A, et al. Prostate specific antigen density: a means of distinguishing benign prostatic hypertrophy and prostate cancer. *J Urol.* 147:815-816, 1992.
- 11 Catalona WJ, Hudson MA, Scardino PT, et al. Selection of optimal prostate specific antigen cut-offs for early detection of prostate cancer: receiver operating characteristic curves. *J Urol.* 152:2037-2042, 1994.
- 12 Smith DS, Catalona WJ. The nature of prostate cancer detected through prostate specific antigen based screening. *J Urol.* 152:1732-1736, 1994.
- 13 Catalona, et al. Evaluation of Percentage of Free Serum Prostate-Specific Antigen to Improve Specificity of Prostate Cancer Screening. *JAMA* 274:1214, 1995.
- 14 Vashi AR, Wojno KJ, Henricks W, et al: Determination of the "reflex range" and appropriate cutpoints for percent free-PSA in 413 men referred for prostatic evaluation using the AxSYM system. *Urology* 49:19-27, 1997.

- 15 Bangma, CH, Kranse R, Blijenberg B, et al: The value of screening tests in the detection of prostate cancer. Part I: Results of a retrospective evaluation of 1726 men. *Urol.* 46:773-778, 1995.
- 16 Yemoto CM, Nolley R, Prestigiacomo AF, et al: Free (f) and total (t) PSA density in patients with prostate cancer (CaP) and benign prostatic hyperplasia (BPH). *J Urol* 155:347A, 1996.
- 17 Bangma, CH, Kranse R, Blijenberg B, et al: The value of screening tests in the detection of prostate cancer. Part II: Retrospective analysis of free-total prostate-specific analysis ration, age-specific reference ranges, and PSA density. *Urol* 46:779-785, 1995.

## TECHNICAL CONSULTATION

Call or Write: **BioCheck, Inc.**  
 323 Vintage Park Dr.  
 Foster City, CA 94404  
 Tel: (650) 573-1968 Fax: (650) 573-1969

080803