

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: T4 Standard/Sample Diluent

Catalog Number: 30004-1 (MB-10102)

Lot Number: RN-

Quantity: _____ml

Product Form: In human serum containing 1.0% (v/v) ProClin-300 as preservative.

Storage: At 2-8°C. **DO NOT FREEZE.**

Suggested Use: For use in the dilution of T4 Standards or patient samples in T4 EIA test.

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.