

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Progesterone-HRP Conjugate

Catalog Number: 70122

Lot Number:

Quantity: _____ ml

Product Form: 1:1 (v/v) (20 mM KPO₄ buffer, 0.85% (w/v) NaCl, pH = 7.2) /glycerol

Storage: -20°C

EIA Titer: Greater than 1:10,000

Suggested Use: For Progesterone Competitive Immunoassay.