

BIOCHECK, INC.
323 Vintage Park Drive
Foster City, CA 94404
Tel: 650-573-1968 Fax: 650-573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Progesterone, Neat Serum

Catalog Number: 70123

Antigen Source: Progesterol-11-Carrier Protein

Lot Number:

Quantity: _____ ml (Neat Serum)

Product Form: Rabbit Anti-Progesterone, Neat Serum, 0.1% NaN₃

Storage: -20°C

Cross-Reactivity:

Progesterone	100%
Androstenedione	0.0%
Androsterone	0.0%
Corticosterone	1.3%
Cortisol	0.0%
Cortisone	0.0%
11-Deoxycortisol	0.0%
11-Deoxycorticosterone	1.5%
Dihydroepiandrosterone	0.0%
Dihydroepiandrosterone sulfate	0.0%
20-Dihydroprogesterone	0.0%
Dihydrotestosterone	0.0%
Epiandrosterone	0.0%

16-Epiestriol	0.0%
Estradiol	0.0%
Estriol	0.0%
Estriol-3-glucuronide	0.0%
Estriol-16-glucuronide	0.0%
Estriol-16-sulfate	0.0%
Estrone	0.0%
5-Pregnandione	2.2%
Pregnenolone	0.4%
17a-OH-Pregnenolone	0.1%
17a-OH-Progesterone	1.9%
Prednisone	0.0%
Testosterone	0.3%

Suggested Use: For Competitive Progesterone Immunoassay.