

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Testosterone-HRP Conjugate

**Catalog Number:** 70124

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** 1:1 (v/v) (20 mM KPO<sub>4</sub> buffer, 0.85% (w/v) NaCl, pH = 7.2) /glycerol

**Storage:** -20°C

**EIA Titer:** Greater than 1:10,000

**Suggested Use:** For Testosterone Competitive Immunoassay.