

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Serum IgG (H+L), DEAE Purified

**Catalog Number:** 70149 (2001)

**Lot Number:** \_\_\_\_\_

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Form:** Liquid

**Purity:** Greater than 98% pure

**Buffer:** 15 mM potassium phosphate buffer, 0.85% sodium chloride, pH=7.4

**Preservative:** 0.10% sodium azide  
This product is filtered through a 0.2 micron membrane.

**Storage:** 2-8°C