

**PRODUCT SPECIFICATIONS**

**Product Name:** Mouse Serum IgG (H+L), DEAE Purified

**Catalog Number:** 70149

**Lot Number:** \_\_\_\_\_

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Form:** Liquid

**Purity:** Greater than 98% pure

**Buffer:** 15 mM potassium phosphate buffer, 0.85% sodium chloride, pH=7.4

**Preservative:** 0.10% sodium azide  
This product is filtered through a 0.2 micron membrane.

**Storage:** 2-8°C