

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-Human CA 19-9  
Tumor Marker (GSLA-3)

**Catalog Number:** 70183 (GSLA-3)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg ( = \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Subclass:** IgG<sub>1</sub>

**Product Form:** Protein A column chromatography purified immunoglobulin fraction in 0.015M KPO<sub>4</sub> buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN<sub>3</sub>.

**Specificity:** Human GI Tumor Antigen (CA19-9)

**Affinity Constant:**  $8 \times 10^9$  l/mol.

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

**Cross-Reactivity:** No cross-reactivity with Human AFP, CEA, PAP, PSA, CA-125, or CA15-3.

**Suggested Use:** For immunohistochemistry (IHC) staining