

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human IgE

Catalog Number: 70188 (121)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column purified immunoglobulin fraction in 0.015 M KPO_4 buffer, 0.85% NaCl, 0.1% NaN_3 , pH = 7.40

Specificity: Human IgE

Cross-Reactivity: No cross reactivity with human IgA, IgD, IgG, IgM

Suggested Use: For plate coating in the human IgE sandwich immunoassay.

Storage: At 2-8°C. DO NOT FREEZE.

Pairing Information: This antibody pairs with goat anti-human IgE-HRP Conjugate (P/N: 80018) in the human IgE sandwich immunoassay.