

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human β -2 Microglobulin

Catalog Number: 70191

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG1

Product Form: Protein A column chromatography purified immunoglobulin fraction in 0.015 M KPO_4 buffer, pH=7.40, containing 0.85% NaCl and 0.1% NaN_3 .

Specificity: Human β -2 Microglobulin

Storage: At 2-8° C, DO NOT FREEZE. Precipitation may occur upon freezing.

Suggested Use: For solid phase coating (capture).

Pairing Information: This antibody pairs with sheep polyclonal Anti-Human β -2 Microglobulin (affinity purified, P/N 70320) in the human β -2 microglobulin sandwich immunoassay.