

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Goat Polyclonal Anti-Human Alpha Fetoprotein (AFP), Affinity Purified
Catalog Number:	70193
Antigen Source:	Human Cord Blood AFP
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Product Form:	AFP-Affinity column chromatography purified immunoglobulin fractions in 0.015M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .
Specificity:	Human Alpha Fetoprotein (AFP)
Storage:	At 2-8° C. <u>DO NOT FREEZE</u> . Precipitation may occur upon freezing.
Cross-Reactivity:	No cross-reactivity with Human Albumin, CEA, Ferritin, hCG, PAP, and PSA.
Suggested Use:	For solid phase immobilization (capture) or enzyme conjugation (detection) in the human AFP sandwich immunoassay. This antibody pairs with AFP MoAb (P/N: 70448).