

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Anti-Human Cardiac Myoglobin, Affinity Purified.

Code: 70196

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: Myoglobin affinity column chromatography purified antibodies in 10 mM sodium phosphate buffer, pH = 7.20, containing 0.85% NaCl and 0.85% NaN_3 .

Storage: 2-8° C

Specificity: This antibody reacts with cardiac and skeletal forms of myoglobin.

Suggested Use: For Enzyme Conjugation (Detection).

Pairing Information: This antibody pairs with myoglobin MoAb (P/N: 70131) in the cardiac myoglobin sandwich immunoassay.