

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human Cardiac Myoglobin,
Affinity Purified.

Code: 70196

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: Myoglobin antigen affinity column chromatography
purified antibodies in 15 mM potassium phosphate
buffer, pH = 7.40, containing 0.85% NaCl and 0.1%
NaN₃.

Storage: 2-8° C

Specificity: This antibody reacts with cardiac and skeletal forms of
myoglobin.

Suggested Use: For enzyme conjugation (detection).

Pairing Information: This antibody pairs with human myoglobin MoAb
(P/N: 70131) in the human cardiac myoglobin sandwich
immunoassay.