

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Anti-CK-MM, Affinity Purified

**Catalog Number:** 70225

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** Liquid

**Purity:** >95% by SDS-PAGE

**Buffer:** Affinity column chromatography purified antibody in 10 mM sodium phosphate buffer, pH= 7.2, containing 0.85% NaCl and 0.05%  $NaN_3$ .

**Preservative:** 0.05%  $NaN_3$

**Storage:** At 2-8°C.

**Suggested Use:** **For Enzyme Conjugation (Detection)**

**Pairing Information:** This antibody pairs with CK-MB MoAb (P/N: 70438) in the CK-MB sandwich immunoassay.