

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Anti-Human CK-MM, Affinity Purified

**Catalog Number:** 70225

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** Liquid

**Purity:** >95% by SDS-PAGE

**Buffer:** CK-MM antigen affinity column chromatography purified antibody in 10 mM sodium phosphate buffer, pH= 7.40, containing 0.85% NaCl and 0.1%  $NaN_3$ .

**Storage:** At 2-8°C.

**Suggested Use:** For enzyme conjugation (detection).

**Pairing Information:** This antibody pairs with CK-MB MoAb (P/N: 70438 or 70444) in the human CK-MB sandwich immunoassay.