

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human C-Reactive Protein (CRP) (Total IgG Fractions)

Catalog Number: 70243

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.015 M KPO_4 buffer, pH=7.40 containing 0.85% NaCl, 0.1% NaN_3 .

Storage: At 2-8° C, DO NOT FREEZE.

Specificity: Human C-Reactive Protein.

Suggested Use: For solid phase coating (capture).

Pairing Information: This antibody pairs with goat polyclonal anti-human CRP (affinity purified, P/N: 70189) or with rabbit polyclonal anti-CRP (affinity purified, P/N: 70648) in the human CRP sandwich immunoassay.