

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Hepatitis B Virus Surface Antigen (HBsAg)

Catalog Number: 70254 (G₁E₁₂)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG₁

Product Form: DEAE column chromatography purified immunoglobulin fraction in 0.015 M KPO₄ buffer, pH=7.40 containing 0.85% NaCl, 0.1% NaN₃

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Reacts with all 4 subtypes, adr, adw, ayr, and ayw of the hepatitis B virus surface antigen.

Affinity Constant: 9×10^9 l/mol.

Suggested Use: For enzyme conjugate or solid phase coating.

Pairing Information: This antibody pairs with Goat Anti-HBsAg Antibody (affinity purified) or with HBsAg MoAb (P/N: 70277) in the HBsAg sandwich immunoassay.