

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Human Chorionic Gonadotropin Antigen  
Iodination Grade

**Catalog Number:** 70264 (4002)

**Lot Number:**

**Source:** Human Pregnancy Urine.

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Specific Immunological Activity:** \_\_\_\_\_ IU/mg/ml

**Reference Preparation:** WHO, 3<sup>rd</sup> IS (75/537)

**Purity:** > 98%. As determined by SDS-PAGE.

**Product Form:** In 0.015M KPO<sub>4</sub> buffer, 0.85% NaCl,  
0.1% NaN<sub>3</sub>, pH=7.40

**Level of Specific Contaminants:** hCG 100%  
hLH < 0.01%  
hFSH < 0.01%  
hTSH < 0.02%

**Storage:** - 20°C or below.