

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Human Alpha Fetoprotein (AFP) (Total IgG Fractions)

Catalog Number: 70269 (# 08)

Antigen Source: Human Cord Blood AFP

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.015M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Specificity: Human Alpha Fetoprotein (AFP)

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Cross-Reactivity: No cross-reactivity with Human Albumin, CEA, Ferritin, hCG, PAP, and PSA.

Suggested Use: For solid phase immobilization in the AFP sandwich immunoassay.

Pairing Information: This antibody pairs with AFP MoAb-HRP Conjugate (P/N: 80005) **in the AFP sandwich immunoassay.**