

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICAITONS

Product Name: Testosterone-HRP Conjugate

Catalog Number: 70287-HRP

Lot Number:

Quantity: _____ ml

Product Form: 1:1 (v/v) (20 mM M KPO₄, 0.85% w/v NaCl, pH = 7.2)
/ glycerol

Storage: -20°C

EIA Titer: Greater than 1:10,000

Suggested Use: For competitive Testosterone immunoassay.