

BIOCHECK, INC.
323 Vintage Park Drive
Foster City, CA 94404
Tel: 650-573-1968 Fax: 650-573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Testosterone, Neat Serum

Catalog Number: 70287

Antigen Source: Testosterone-Carrier Protein

Lot Number:

Quantity: _____ ml (Neat Serum)

Product Form: Rabbit Anti-Testosterone, Neat Serum, 0.1% NaN₃

Storage: -20°C

Cross-Reactivity:

Testosterone	100%	Estriol	0.0%
Androstenedione	0.0%	Estriol-3-glucuronide	0.0%
Androsterone	0.0%	Estriol-16-glucuronide	0.0%
Dihydroepiandrosterone	0.0%	Estriol-16-sulfate	0.0%
Dihydrotestosterone	5.0%	Estrone	0.0%
Epiandrosterone	0.0%	17a-Pregnenolone	0.0%
16-Epiestrol	0.0%	17a-Progesterone	0.0%
Estradiol	0.0%	Progesterone	0.0%

Suggested Use: For Competitive Testosterone Immunoassay.