

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Polyclonal Anti-Human Alpha Fetoprotein (AFP)  
(Total IgG Fractions)

**Catalog Number:** 70295 (NS-AFP)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg ( = \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** DEAE column chromatography purified immunoglobulin fractions in 0.015M  $KPO_4$  buffer, pH = 7.40, containing 0.85% NaCl, 0.1%  $NaN_3$ .

**Specificity:** Human Alpha Fetoprotein (AFP)

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

**Cross-Reactivity:** No cross-reactivity with Human Albumin, CEA, Ferritin, hCG, PAP, and PSA.

**Suggested Use:** For solid phase immobilization in the AFP sandwich immunoassay

**Pair Information:** This antibody pairs with AFP MoAb-HRP Conjugate (P/N: 80005) **in the AFP sandwich immunoassay.**