

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human GI Tumor Marker (CA 19-9)

Catalog Number: 70313

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG₁

Product Form: Protein A column chromatography purified immunoglobulin fraction in 0.015M KPO₄ buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN₃. Ab solution filtered through a 0.2 μ filter.

Specificity: Human GI Tumor Antigen (CA19-9)

Affinity Constant: N/A

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Cross-Reactivity: No cross-reactivity with Human AFP, CEA, PAP, PSA, CA-125, or CA15-3.

Suggested Use: **For enzyme conjugation and solid phase Coating. This same antibody is for both purpose.**

Pairing Information: This antibody pairs with itself in the CA 19-9 sandwich immunoassay.