

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human Prolactin
Catalog Number:	70317 (172)
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Product Form:	DEAE column chromatography purified immunoglobulin fractions in 0.015M KPO ₄ buffer, pH=7.40, containing 0.85% NaCl, 0.1% NaN ₃ .
Specificity:	Human prolactin
Cross-Reactivity:	No cross reactivity with hCG, LH, FSH, growth hormone and TSH.
Affinity Constant:	1 x 10 ¹⁰ l/mol.
Storage:	At 2-8° C. Do Not Freeze. Precipitation may occur upon freezing.
Suggested Use:	For enzyme conjugation in the human Prolactin sandwich immunoassay.
Pairing Information:	This antibody pairs with Prolactin MoAb (P/N: 70336) in <u>the Prolactin sandwich immunoassay.</u>