

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Sheep Polyclonal Anti-Human β -2 Microglobulin, Affinity Purified

Catalog No: 70320

Lot Number:

Immunogen: Highly purified human β -2 Microglobulin

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: Affinity column chromatography purified sheep polyclonal anti-human β -2 microglobulin in 10 mM sodium phosphate buffer, pH = 7.20, containing 0.85% NaCl and 0.05% NaN₃.

Storage: 2-8° C

Specificity: This antibody reacts with human β -2 Microglobulin.

Suggested Use: For Enzyme Conjugation (Detection)

Pairing Information: This antibody pairs with β -2 microglobulin MoAb (70191) in the **human β -2 microglobulin** sandwich immunoassay.