

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human CEA

Catalog Number: 70321

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: Protein A column chromatography purified immunoglobulin fractions in 0.015M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Specificity: Human CEA

Affinity Constant: 3×10^{10} l/mol.

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Cross-Reactivity: No cross-reactivity with human NCA, NCA2 or BGP1 (Biliary glycoprotein-I)

Suggested Use: For immunohistochemistry (IHC) staining