

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	<b>Goat Polyclonal Anti-Human CEA, Affinity Purified</b>
<b>Catalog Number:</b>	70322
<b>Antigen Source:</b>	Metastasized human liver tumors
<b>Lot Number:</b>	
<b>Quantity:</b>	_____ mg (= _____ ml)
<b>Antibody Concentration (<math>A_{280nm}</math>):</b>	_____ mg/ml
<b>Product Form:</b>	CEA antigen affinity column chromatography purified immunoglobulin fractions in 0.015 M $KPO_4$ buffer, pH=7.40 containing 0.85% NaCl, 0.1% $NaN_3$
<b>Storage:</b>	At 2-8° C
<b>Specificity:</b>	Human CEA
<b>Cross-Reactivity:</b>	No cross-reactivity with human NCA, NCA2 or BGP1 (Biliary glycoprotein-I).
<b>Suggested Use:</b>	For enzyme conjugation (detection) in the CEA sandwich immunoassay
<b>Pairing Information:</b>	This antibody pairs with mouse monoclonal anti-human CEA (P/N: 70458) in the human CEA sandwich immunoassay.