

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human CEA, Affinity Purified

Catalog Number: 70322 (3005/379)

Antigen Source: Metastasized human liver tumors

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: CEA Affinity column chromatography purified immunoglobulin fractions in 0.015 M KPO_4 buffer, pH=7.40 containing 0.85% NaCl, 0.1% NaN_3

Storage: At 2-8° C

Specificity: Human CEA

Cross-Reactivity: No cross-reactivity with human NCA, NCA2 or BGP1 (Biliary glycoprotein-I).

Suggested Use: **For enzyme conjugation in the CEA sandwich immunoassay**

Pairing Information: This antibody pairs with mouse monoclonal anti-CEA (P/N: 70357 and 70370) in the CEA sandwich immunoassay.