

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti Ferritin,
total IgG Fractions, DEAE Purified

Catalog Number: 70326 (R-9293)

Antigen Source: Human liver and spleen ferritin

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin
Fractions in 0.015 M KPO_4 buffer, pH=7.40 containing
0.85% NaCl, 0.1% NaN_3 .

Storage: At 2-8°C.

Specificity: Human liver and spleen ferritin

Cross-Reactivity: No cross-reactivity with human albumin, AFP, CEA,
hCG, PAP and PSA.

Suggested Use: For solid phase immobilization in the ferritin sandwich
immunoassay.

Pairing Information: This antibody pairs with ferritin MoAb (P/N: 70166) in the
ferritin sandwich immunoassay.