

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Sheep Polyclonal Anti-Human Growth Hormone, total IgG, DEAE Purified

Catalog Number: 70332 (PS-602)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.015 M KPO_4 buffer, pH = 7.40 containing 0.85% NaCl, 0.1% NaN_3 .

Storage: 2-8°C. DO NOT FREEZE.

Specificity: Human Growth Hormone

Cross-Reactivity: No cross-reactivity with hCG, hLH, hFSH, hProlactin or hTSH.

Suggest Use: For solid phase coating or for further affinity purification.

Pairing Information: This antibody pairs with HGH MoAb (P/N: 70386) in the sandwich immunoassays for human growth hormone.