

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Rabbit Polyclonal Anti-Human Prostate Specific Antigen Total PSA (Total IgG Fractions)

**Catalog Number:** 70339 (1110)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280}$  nm):** \_\_\_\_\_ mg/ml

**Product Form:** DEAE column chromatography purified immunoglobulin fraction in 0.015M  $KPO_4$  buffer, pH = 7.40, containing 0.85% NaCl, 0.1%  $NaN_3$ .

**Specificity:** Human Prostate-specific antigen

**Cross-Reactivity:** No cross-reactivity with Human Albumin, AFP, CEA, Ferritin, hCG, PAP, CA 125, CA 15-3, or CA 19-9.

**Storage:** At 2-8° C.

**Suggested Use:** For solid phase coating.

**Pairing Information:** This antibody pairs with mouse monoclonal anti-PSA (P/N: 70338) in the Total PSA sandwich immunoassay.