

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human-Beta-TSH,
Affinity Purified

Catalog Number: 70342 (G-7003)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: TSH and alpha hCG affinity columns purified immunoglobulins in 0.015 M KPO_4 buffer, pH=7.40, containing 0.85% NaCl and 0.1% NaN_3 .

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Human Beta-TSH

Cross-Reactivity: Less than 0.1% cross-reactivity with alpha hCG, LH, FSH, growth hormone and prolactin.

Suggested Use: For enzyme conjugation or solid phase coating.

Pairing Information: This antibody pairs with TSH MoAb (P/N: 70204)
in the TSH sandwich immunoassay.