

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Alpha-TSH, Affinity Purified

Catalog Number: 70343 (G-7005)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: Alpha hCG affinity columns purified Immunoglobulins from goat anti-TSH serum in 0.015 M KPO_4 buffer, pH=7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Purity: Greater than 98% by protein gel electrophoresis

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Human Alpha-TSH

Cross-Reactivity: Cross react with Alpha hCG, Alpha-LH, and Alpha-FSH.

Suggested Use: For enzyme conjugation or solid phase coating