

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Polyclonal Anti-Rabbit IgG (H&L),  
Affinity Purified

**Catalog Number:** 70347 (3002)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg ( = \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Form:** Rabbit IgG affinity column purified goat anti-rabbit  
IgG(H+L). In 0.10 M potassium phosphate buffer,  
pH = 7.80 containing 0.05%  $NaN_3$ .

**Purity:** Greater than 98% pure

**Storage:** 2-8°C, DO NOT FREEZE

**Specificity:** This antibody reacts with heavy chains on rabbit IgG, and  
with light chains on all rabbit immunoglobulins. No  
reactivity is observed to non-immunoglobulin rabbit  
serum proteins.

**Note:** Slight precipitation of antibody may occur on storage at  
2-8°C. Centrifuge or refilter if necessary.