

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Rabbit IgG (H+L),
Affinity Purified

Catalog Number: 70347 (3002)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Form: Rabbit IgG affinity column purified goat polyclonal anti-rabbit IgG (H+L). In 0.10 M potassium phosphate buffer, pH = 7.80 containing 0.05% NaN_3 .

Purity: Greater than 98% pure

Storage: 2-8°C, DO NOT FREEZE

Specificity: This antibody reacts with heavy chains on rabbit IgG, and with light chains on all rabbit immunoglobulins. No reactivity is observed to non-immunoglobulin rabbit serum proteins.

Note: Slight precipitation of antibody may occur on storage at 2-8°C. Centrifuge or refilter if necessary.