

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Anti-Fatty Acid Binding Protein (FABP),
Affinity Purified

Catalog Number: 70385

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: Affinity purified column chromatography purified immunoglobulin fraction in 0.015 M KPO_4 , pH = 7.40, containing 0.85% NaCl and 0.1% NaN_3 .

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Human Fatty Acid Binding Protein

Suggested Use: For solid phase coating and HRP-conjugation.