

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human Growth Hormone

Catalog Number: 70386 (A36020047A)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG₁

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.015M KPO₄ buffer, pH = 7.40, containing 0.85% NaCl and 0.1% NaN₃.

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Human Growth Hormone

Affinity Constant: 8×10^9 l/mol.

Suggested Use: For enzyme conjugation in the **human Growth Hormone sandwich immunoassay**.

Pairing Information: This antibody pairs with Sheep Anti-Human Growth Hormone (P/N: 70332) in the Human Growth Hormone Sandwich Immunoassay.