

BIOCHECK, INC.
323 Vintage Park Drive
Foster City, CA 94404
Tel: 650-573-1968 Fax: 650-573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Estradiol, Neat Serum

Catalog Number: 70388

Antigen Source: Estradiol-6-CMO-BSA

Lot Number:

Quantity: _____ ml (Neat Serum)

Product Form: Rabbit Anti-Estradiol, Neat Serum, 0.1% NaN₃.

Storage: -20°C or below.

Cross-Reactivity:

Estradiol-17	100%	Estradiol-3-glucuronide	0.00%
Androsterone	0.00%	Estriol	0.05%
Androstenedione	0.00%	Estriol-16-glucuronide	0.00%
Corticosterone	0.00%	Estrone	0.20%
Cortisone	0.00%	Estrone-3-sulfate	0.00%
11-Deoxycortisol	0.00%	11-Hydroxyprogesterone	0.00%
11-Deoxycorticosterone	0.00%	17a-Hydroxyprogesterone	0.00%
21-Deoxycortisol	0.00%	17a-Pregnenolone	0.00%
Dihydroepiandrosterone	0.00%	17a-Progesterone	0.00%
Dihydrotestosterone	0.00%	Pregnanediol	0.00%
20-Dihydroprogesterone	0.00%	Pregnanetriol	0.00%
Epiandrosterone	0.00%	Pregnenolone	0.00%
16-Epiestriol	0.00%	Progesterone	0.00%
Estradiol-3-sulfate	0.00%	Testosterone	0.00%

Suggested Use: For Estradiol Competitive Immunoassay.