

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Goat Polyclonal Anti-Whole hCG, Affinity Purified
Catalog Number:	70396
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Product Form:	In 0.1 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN_3 , pH = 7.80
Purity:	Greater than 98% pure by PAGE
Preservative:	0.05% sodium azide, filtered through a 0.2 micron membrane.
Storage:	At 2-8°C.
Specificity:	This antibody reacts with alpha chain and beta chain of the whole hCG molecule.
Pairing Information:	This antibody pairs with free β -hCG MoAb (#70205) in the <u>free β-hCG sandwich immunoassay</u> .
Note:	Slight precipitation of antibody may occur on storage at 2-8°C. If necessary, refilter or centrifuge at 3,500 rpm for 15 minutes to remove the tiny precipitate.