

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human Ovarian Tumor Marker (CA 125)

Catalog Number: 70412

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG₁

Product Form: Protein A column chromatography purified from mouse ascites fluid, provided in PBS buffer, pH 7.5 with 0.05% NaN₃

Specificity: CA125: 100%, Epitope specificity "M-11 like" Group B

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Suggested Use: For enzyme conjugation (detection) use

Pairing Information: This antibody pairs with CA 125 MoAb (P/N: 70400) in the human CA 125 sandwich immunoassay.