

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Human Prostate Specific Antigen (PSA)

Catalog Number: 70413 (4001)

Lot Number:

Source: Human Seminal Fluid

PSA Concentration: _____ mg/ml

Quantity: _____ mg (= _____ ml)

Product Form: Single homogenous batch, sterile filtered through a 0.2 μ membrane in a 0.05 M Phosphate buffer containing 0.15M NaCl and 0.1% NaN₃. pH=7.50

Purity: > 95% (by SDS-PAGE)

Storage: At 2-8°C. DO NOT FREEZE.