

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Goat Polyclonal Anti-Human Prostate Specific Antigen (PSA), Affinity Purified
Catalog Number:	70415
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Form:	PSA antigen-affinity column chromatography purified immunoglobulin fraction in 0.015 M KPO ₄ buffer, pH = 7.40 containing 0.85% NaCl, 0.1% NaN ₃ .
Purity:	> 95% by SDS-PAGE
Storage:	At 2-8°C. <u>DO NOT FREEZE.</u>
Specificity:	Human prostate-specific antigen
Cross-Reactivity:	No cross-reactivity with human albumin, AFP, CEA, Ferritin, hCG, PAP, CA 125, CA 19-9, or CA 15-3.
Suggested Use:	For enzyme conjugation (detection).
Pairing Information:	<ol style="list-style-type: none">1. This antibody pairs with mouse monoclonal anti-PSA (P/N: 70338) for Total PSA sandwich immunoassay.2. This antibody also pairs with mouse monoclonal anti-Free PSA (P/N: 70340) for Free PSA sandwich immunoassay.