

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-PSA, Affinity Purified

Catalog Number: 70423 (3029)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: PSA Affinity column chromatography purified immunoglobulin fraction in 0.015M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Specificity: Human Prostate-specific antigen

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Cross-Reactivity: No cross-reactivity with human albumin, AFP, CEA, Ferritin, hCG, PAP, CA125, CA15-3, or CA19-9.