

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-Hepatitis B Virus e Antigen (HBeAg)

**Catalog Number:** 70427 (C4)

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Subclass:** IgG<sub>1</sub>

**Product Form:** Protein A column chromatography purified immunoglobulin fraction in 0.015 M KPO<sub>4</sub> buffer, pH=7.40 containing 0.85% NaCl, 0.1% NaN<sub>3</sub>

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

**Specificity:** Reacts with the hepatitis B virus e antigen

**Affinity Constant:**  $8 \times 10^9$  l/mol.

**Suggested Use:** For enzyme conjugate, radio labeling, or solid phase coating.

**Pairing Information:** This antibody pairs with HBe MoAb (P/N: 70426) in the HBeAg ELISA sandwich immunoassay.