

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-hCG-Alpha, Affinity Purified

Catalog Number: 70435

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.1 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN₃, pH = 7.80

Purity: Greater than 98% pure by PAGE

Storage: At 2-8°C.

Specificity: This antibody reacts with alpha chain of hCG molecule.

Pairing Information: This antibody pairs with hCG MoAb (P/N: 70330) in the **whole hCG sandwich immunoassay**.

Note: Slight precipitation of antibody may occur on storage at 2-8°C. If necessary, refilter or centrifuge at 3,500 rpm for 15 minutes to remove the tiny precipitate.