

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Chloramphenicol (CAP),
DEAE Purified (Total IgG Fractions)

Catalog Number: 70442

Antigen Source: Chloramphenicol-KLH

Lot Number:

Amount: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.015 M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Specificity: Chloramphenicol (CAP)

Storage: Short Term: 2-8°C (one month)
Long Term: -20°C or below.

Suggested Use: For chloramphenicol Immunoassay