

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Clenbuterol-HRP Conjugate

**Catalog Number:** 70443-HRP

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Product Form:** 1:1 (v/v) (20 mM M KPO<sub>4</sub>, 0.85% w/v NaCl, pH = 7.2) / glycerol

**Storage:** -20°C

**Suggested Use:** For competitive clenbuterol immunoassay