

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Thyroxin (T4)

Immunogen: T4-BSA

Catalog Number: 70446

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.05 M PBS, pH = 7.40, containing 0.1% NaN_3 .

Affinity Constant: 1×10^9 l/mol.

Subclass: IgG₁

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Cross-Reactivity: T2: < 1%
T3: < 1%
rT3: < 1%