

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-Myeloperoxidase (MPO) Antibody

**Catalog Number:** 70457

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Subclass:** IgG<sub>1</sub>

**Product Form:** Protein A column chromatography purified immunoglobulin fraction in 0.015M KPO<sub>4</sub>, pH = 7.40, containing 0.85% NaCl and 0.1% NaN<sub>3</sub>.

**Storage:** At 2-8°C, Do Not Freeze. Precipitation may occur upon freezing.

**Specificity:** Human Myeloperoxidase.

**Suggested Use:** For solid phase coating.

**Pairing Information:** This antibody pairs with mouse monoclonal anti-MPO (P/N:70356) in the **myeloperoxidase** sandwich immunoassay.