

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Human Monoclonal Anti-NT-proBNP

**Catalog Number:** 70464 (2-D)

**Lot Number:** \_\_\_\_\_

**Quantity:** \_\_\_\_\_ mg

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** In 0.015M  $KPO_4$  buffer, pH = 7.40, containing 0.85% NaCl, 0.1%  $NaN_3$

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.